

Technical Data Sheet

Product Description – 90 gram White glassine, silicone C1S

Test Item	Test Method	Unit	Target	Tolerance
Basis Weight	GLI 11	gsm	90	±5
Caliper	GLI 12	micras	76	±3
Release Force* C1S	CLI 13	g/inch	5-20	±6
Tensile MD	ISO 1924-2	kN/m	10	±0.2

Features:

- Solvent-less Platinum catalyzed silicone with easy release silicone on one side.
- Glassine paper provides a smooth surface with excellent caliper control.
- Paper can be direct coated with adhesive.
- Primarily used for labels.

REMARKS

² Typical values will not to be perceived as specification and is only for reference.

I Once the liner is taken into production, user is responsible for determining whether it is fit for a particular purpose, method of application and service life of product.

² This release liner is specially designed and intended for single use purpose only. MexRelease will not guarantee the liner performance on multiple usages.

*Test method FTM10 23°C, TESA 7475 tape, pulled 300mm/minute at 180 degrees



Azucena no. 25, Colonia Ex ejido San Lucas Patoni Tlalnepantla de Baz, Estado de México, C.P 54100

RFC:CME-210317-UA2